
Pellin-Broca-Prismen

Pellin-Broca Prisms

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These prisms are solid dispersing prisms with 90° deviation. They can be envisaged as a combination of two 30° dispersing prisms with one right angle reflection prism. The combination of the two 30° prisms works like a 60° dispersing prism and the 90° prism only does a pure reflection without dispersion.

The two 30° prisms (i.e. the 60° dispersing prism) are transmitted in the so-called „minimum of deviation“ by light of one wavelength. The entire prism deflects light of this wavelength by exactly 90° ; for this wavelength the spectroscopic angular dispersion is maximum.

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If the prism is rotated about axis A (see sketch), then a different wavelength is deviated by 90°. This way the prism can be tuned without the necessity to alter anything else in the geometry of the experimental setup.

The three polished surfaces of the prism produced from one block have a surface accuracy of $\lambda/10$. The angular errors are about 2'.

The following table lists a selection of Pellin-Broca prisms including individual beam cross-sections, material and reflection value R_s and R_p (per surface) of light polarized parallel or perpendicular to the plane of incidence (at some wavelengths λ).

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Other sizes and materials on request

Literatur:

W. M. McClain, Appl. Opt. **12**, (1973), S. 153ff „How to Mount a Pellin-Broca Prism for Laser Work“

J.Reader, Appl. Opt. **12**, (1973), S. 1405, „Comment on: How to Mount a Pellin-Broca Prism for Laser Work“

Literature:

W. M. McClain, Appl. Opt. **12**, (1973), p. 153ff, „How to Mount a Pellin-Broca Prism for Laser Work“

J.Reader, Appl. Opt. **12**, (1973), p. 1405, „Comment on: How to Mount a Pellin-Broca Prism for Laser Work“