
Soleil-Babinet-Kompensatoren

Soleil-Babinet Compensators

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The compensators are used to measure states of polarization or to produce certain retardation values: e.g., $\lambda/2$ or $\lambda/4$ for different wavelengths. The clear aperture is 12 mm in diameter. The compensators are supplied with a cylindrical attachment 25 mm in diameter for mounting in graduated rotating mounts or other holders.

The compensators consist of a long wedge $39 \times 13 \text{ mm}^2$ in size and a short wedge ($13 \times 13 \text{ mm}^2$) with identical wedge angle. The thickness of the $13 \times 13 \text{ mm}^2$ compensating plate is matched to provide zero path difference approximately in the center of the long wedge (the zero position). The wedge in a precision slade mount is adjustable by means of a micrometer screw. The adjustment range is about 13 mm on either side of the zero position.

Based on this adjustment range a guideline for the order numbers for the different compensators is presented in the table. To provide exact measurements, the compensators must be calibrated for the respective wavelength.

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
The flatness of the optical faces is about $\lambda/20$ for the compensators made of quartz and $\lambda/10$ for the version in MgF₂.

Directions for adjustment and operation, including the procedure for calibration, are also supplied.

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Quality of material used:

The field of view is almost free of any imperfections in the crystal structure. per the above table

The field of view contains imperfections. RKQ 0

Available on special order:

Compensators can be supplied with calibration for certain wavelengths in the visible spectral range (633 nm, 546 nm, 435 nm).

Calibration per wavelength RKQ 001